

## PINWORMS INVESTIGATION (adhesive spatula)

### INSTRUCTIONS

- The requested examination must be performed before getting up, without mobilizing the intestine or sanitizing.
- The spatula must be applied in the anal margins, using the rubberized side of it, holding it with the hand by the lid, and pressing it several times on the skin around the anal orifice, without introducing it in the anus.
- THIS PROCEDURE MUST BE REPEATED FOR THREE DAYS IN A ROW USING THE SAME SPATULA, SAVING IT INSIDE THE TUBE IN BETWEEN TAKES.
- Label the tube with the patient's full name and member number and deliver it to the Clinical Laboratory from 7:30 a.m. to 6:00 p.m.